

Zimmer PlumeVac Smoke Evacuator - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Zimmer PlumeVac Smoke Evacuator provides flexible smoke evacuation for surgical procedures. The PlumeGuard 35Hr Filter delivers up to 35 hours of performance and features a covered 3-port design that accommodates tubing sizes of 1/4" (6.4mm), 3/8" (9.5mm), and 7/8" (22mm), enabling versatility across different procedures and accessory configurations.

The unit operates quietly and effectively while capturing particles with excellent efficiency. It includes three fan speeds and a filter life meter to monitor performance. The system works with most ESU units and is designed for low cost per procedure.

Package includes the evacuation unit, power cord, user's guide, pneumatic footswitch, and a new PlumeSafe ULPA filter ready for immediate use.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)