

# Alaris Medley 8015 PC - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Alaris Medley 8015 PC Infusion Pump Controller : Refurbished. is the heart of the Alaris® System and can be used for neonatal, pediatric, and adult patients. It allows for up to 4 detachable infusion or monitoring modules. The Guardrails® Software helps minimize dangerous IV medication errors. Additionally, the software tracks and measures system performance as well as acts as an aid in assisting national safety standards.

The 8015 is extremely flexible and can come in a variety of configurations. The 8015 is the controller, the “brain” of the system. It must be accompanied by additional components to administer the therapy needed.

The Alaris Medley 8015 PC can support up to 4 attached modules.

## FEATURES

- For Neonatal, pediatric, and adult patients alike
- Minimize dangerous IV medication errors
- Tracks and measures system performance
- Helps Facility with national safety standards
- Guardrails® Suite MX is the most comprehensive safety software available

## SPECIFICATIONS

- Approximate Dimensions: 6.9"W x 8.8"H x 9" D
- Approximate Weight: 7.2 lbs.
- Battery Operation: Battery run time is a function of the number of modules attached
- Power Requirements: 100-240V, 50/60 Hz, 150VA Max

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)