

Alaris Medley CareFusion 8100 Infusion Pump Module – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Alaris Medley CareFusion 8100 Infusion Pump Module: Refurbished. This Infusion pump provides excellent patient care and is designed for use with the following: Neonatal Pediatrics Adults It is lightweight and can help protect any type of infusion including intermittent and continuous deliveries.

This module can be pole-mounted when coupled with the Alaris 8015 PC and applied in a wide variety of clinical settings. The Alaris 8100 Infusion Pump Module has been Completely Refurbished and Biomedically Inspected.

Note: The Alaris 8100 does not work independently, and must be used in conjunction with Alaris 8000 or 8015 POC units. Up to 4 pumps can be coupled with one controller.

FEATURES

- Large Volume Pumping Module
- Easily Attached to Pole
- Light Weight
- For use with Neonatal, Pediatric, and Adults
- For Delivery of fluids, blood products, and medications

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)