

Baxter Sigma Spectrum Infusion Pump – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Baxter Sigma Spectrum Infusion Pump – Refurbished. The Baxter Sigma Spectrum is a compact, advanced infusion system designed to improve patient safety and clinical efficiency. This refurbished unit delivers the same reliable performance as a new device at a significantly lower cost.

Ideal for hospitals, outpatient centers, and long-term care facilities, this infusion pump meets high standards for safety, usability, and interoperability. Designed with smart infusion technology, the Sigma Spectrum helps reduce medication errors through integrated drug libraries and safety alerts. Its intuitive user interface features a bright, easy-to-read color screen that displays onboard help menus.

This makes it easier for clinicians to operate and troubleshoot in real time.

FEATURES

- Intuitive Design
- Bright LCD Display
- Multiple Safety System
- Standard IV Admin Sets
- Built-in Drug Safety Library

SPECIFICATIONS

- Weight: Approx 2 lbs (depending on accessories)
- Dimensions: 64.2”D x 4.2”W x 6.3”H (excluding pole clamp)
- Volume Range: 0.1 to 99.9 mL (in 0.1 ML increments)
- Volume Range: 100 to 999 ML (in 1.0 mL increments)
- Mounting Clamp (pole)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)