

Braun Vista Basic Infusion Pump – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Braun Vista Basic Infusion Pump – Refurbished. This infusion pump is completely refurbished and bio-medically inspected by Venture Medical. Furthermore, the Braun Vista Basic Infusion Pump provides rapid, intuitive infusions in a small package. This unit uses administration sets for both gravity and pump infusions.

Additionally, smart technology helps prevent pump-related medical errors. Straight-line tubing loads quickly into the pump. Simply enter the rate and volume and press Start. Helpful prompts walk the user through the programming process Easy to use, teach, and learn.

Perfect for homecare settings Pump-based free-flow protection prevents accidental boluses when the door opens. Adjustable occlusion settings help protect patients from vein injuries.

FEATURES

- User Friendly Design
- Bright LCD Display
- Multiple Safety Systems
- Rate Range: 0.1-800 ml/hr in 0.1 ml increments with Vista™ tubing
- Occlusion Pressure: Variable Settings: Low – 8 PSI, High – 17 PSI

SPECIFICATIONS

- Weight: Approx 6.8 lbs
- Dimensions: 5.5"W x 9.4"H x 7.8"D(excluding pole clamp)
- Rate Range: 0.1-800 ml/hr in 0.1 ml increments with Vista™ tubing
- Occlusion Pressure: Variable Settings: Low – 8 PSI, High – 17 PSI
- Volume Range: 0.1-9999.9 ml

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)