

Alaris Medley 8015 + 8100 Combo Infusion Pump System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Alaris Medley 8015 PC + 8100 Combo Infusion Pump System: Refurbished. The most popular combination of pump + controller. In addition, the Alaris® Medley System can be used for neonatal, pediatric, and adult patients.

Furthermore, it allows for up to 4 detachable infusion or monitoring modules. Safety is key, and the Guardrails® Software helps minimize dangerous IV medication errors. Additionally, the software tracks and measures system performance as well as acts as an aid in assisting national safety standards. Furthermore, this IV pump combo is extremely flexible and can come in a variety of configurations.

The 8015 is the controller, the “brain” of the system.

FEATURES

- For Neonatal, pediatric, and adult patients alike
- Minimize dangerous IV medication errors
- Tracks and measures system performance
- Helps Facility with national safety standards
- Guardrails® Suite MX is the most comprehensive safety software available

SPECIFICATIONS

- Approximate Dimensions: 6.9"W x 8.8"H x 9" D
- Approximate Weight: 7.2 lbs.
- Battery Operation: Battery run time is a function of the number of modules attached
- Power Requirements: 100-240V, 50/60 Hz, 150VA Max

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)