

# Hospira ICU Medical Plum 360 Infusion Pump - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Hospira ICU Medical Plum 360 Infusion Pump - Refurbished. The Plum 360 IV pump delivers reliable performance that helps reduce medication errors and improve patient safety. Hospitals and clinics trust its consistent accuracy.

Moreover, it offers industry-leading cybersecurity features to protect patient data. With full IV-EHR interoperability, the Plum 360 seamlessly integrates into modern healthcare systems, enhancing workflow and care coordination. In addition, the Plum 360 uses the innovative PlumSet cassette for automated back-priming.

This feature allows clinicians to remove air from the IV line easily. They can do this without disconnecting from the patient, which maintains a closed system. As a result, infection risks decrease, and setup becomes more efficient.

## FEATURES

- The Plum 360 utilizes the unique PlumSet cassette that features air management that doesn't require disconnecting.
- Use the pump to automatically remove air through the cassette's secondary administration port without disconnecting the line
- Capture 1 mL of air in the cassette's air trap - an industry exclusive feature - before the alarm is activated
- Help eliminate the potential for contaminants by keeping the line connected
- Hospira's Plum infusion systems are unique in the way they connect the secondary line directly to the PlumSet cassette, so the pump can differentiate between primary and secondary infusion lines and deliver each as intended

## SPECIFICATIONS

- Battery Life: Up to 3 hours of use per charge
- Weight: Approx. 10 lbs (4.5 kg)
- Dimensions: 8" H x 8" W x 6" D
- Volume Range: 0.1-9999 mL
- Delivery Rate Accuracy: ±5% (standard range)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*