


MCM Siderail Socket - Aluminum | MCM 161



Stock Photo - Options May Vary

FACTORY NEW

The MCM 161 is an aluminum siderail socket designed for demanding clinical environments. Constructed from durable aluminum, it offers a lighter yet stronger alternative to traditional clamps, improving handling while maintaining a secure hold during surgical procedures and clinical use.

The MCM 161 accepts flat bar accessories, making it highly versatile across different setups. It attaches easily to side rail ends or notches, allowing flexible positioning to meet the specific needs of various medical applications.

Its lightweight construction enhances mobility and efficiency, reducing strain on medical staff during equipment setup. Despite its minimal weight, the MCM 161 delivers dependable strength for repeated use in operating rooms and surgical centers.

FEATURES

- Accepts standard flat bar accessories
- Aluminum construction provides a stronger, light clamp
- Attaches at rail ends or notches

SPECIFICATIONS

- Accepts Standard Flat Bar Accessories that measure 1" x 1/4"
- Construction: Aluminum

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)