

Amsco 3000 Series Re-manufactured 16" Pre-Vac Sterilizer | 3013



Stock Photo - Options May Vary

FACTORY NEW

The Amsco 3000 Series: Re-manufactured 16" Pre-Vac Sterilizer Model Number: 3013. Features include Eagle 3000 Stage 3 Series controls, a single door, and cabinet mounting.

It also features an advanced microcomputer control system that provides the latest standards for cycle setup, selection, and monitoring. Once you set the settings, the microcomputers accurately monitor and control system operations.

Additional features: The sterilizer comes equipped with a new Reimers steam generator, one (1) interior rack, shelving, and a warranty. Note: To view additional models, click the following link, Sterilizers, Disinfectors, & Washers.

FEATURES

- Eagle 3000 Stage 3 Series Controls
- Advanced Microcomputer Control System
- New Reimers Steam Generator
- Accurately Monitors and Controls System Operations
- Single Door

SPECIFICATIONS

- Chamber Dimensions: 16" x 16" x 20"
- Overall Exterior Dimensions (including door lock-wheel): 74.5"H x 31.5"W x 38.5"D
- Controls: Eagle 300 Stage 3 Series
- Steam Generator: New Reimers
- Operating Temperature: 250° - 270°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)