

Amsco 3000 Series Re-manufactured 20" Pre-Vac Sterilizer | 3023



Stock Photo - Options May Vary

FACTORY NEW

The Amsco 3000 Series Re-manufactured 20" Pre-Vac Sterilizer Model Number: 3023. Features include Eagle 3000 Stage 3 Series controls, a single door, and cabinet mounting.

It also features an advanced microcomputer control system that provides the latest standards for cycle setup, selection, and monitoring. Once you set the settings, the microcomputers accurately monitor and control system operations.

Additional features: The sterilizer includes a new Reimers steam generator, an interior rack and shelving assembly, and a warranty. Note: To view additional models, click the following link, Sterilizers, Disinfectors, & Washers.

FEATURES

- Eagle 3000 Stage 3 Series controls
- Single Door
- Cabinet Mounted
- Operating Temperature of 250 - 270°F
- Advanced Microcomputer Control System

SPECIFICATIONS

- Chamber Dimensions: 20" x 20" x 38"
- Overall Exterior Dimensions (including door lock-wheel): 74.5"H x 35.5"W x 50.5"D
- Operating Temperature: 250 - 270°F
- Exterior Construction: Stainless Steel
- Warranty: 1 Year Parts and Labor, 8 Year Vessel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)