

# Amsco/Steris Washer/Disinfector, Re-manufactured | 444



Stock Photo - Options May Vary

## FACTORY NEW

The Re-manufactured Amsco/Steris Washer/Disinfector, Model Number: 444. Features Include Double Air Actuated Power Doors, Stage 3000 Controls, Electric Steam Configuration (480v 3 Phase), or House Steam and Pure Water Rinse. It also features a Non-Vented Drying System and a Thermal Printer for Accurate Record Keeping.

Additional Features: This Unit comes with (1) Manifold Loading Rack, (1) Transfer Cart Workstation, and also a 1 Year Parts and Labor Warranty. The team completely disassembles the unit to the bare core, rebuilds it using all OEM parts and procedures, and conducts a biomedical inspection.

Note: To view additional models, click the following link, Sterilizers, Disinfectors, & Washers.

## FEATURES

- Double Air Actuated Power Doors
- Stage 3000 Controls: Electric Steam Configuration or House Steam Configuration
- Pure Water Rinse
- Non-Vented Drying System
- Thermal Printer for Accurate Record Keeping

## SPECIFICATIONS

- Overall Exterior Dimensions: 42"W x 80"H x 32 5/8"D
- Interior Chamber Load Capacity: 24" x 26" x 24"
- Electric Steam Configuration: 480v 3 Phase
- House Steam Configuration: 208v or 480v
- Exterior Cabinet: Stainless Steel

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)