

Blickman Built-In Warming Cabinet | SW30TG



Stock Photo - Options May Vary

FACTORY NEW

Blickman Built-In Warming Cabinet SW30TG delivers precise heating performance for modern surgical environments. This unit fits recessed or built-in installations, and it works perfectly in operating room wall consoles. It warms blankets or solutions with controlled temperatures from 90°F to 160°F.

The vibrant graphical display uses 160 x 180 pixels in dual colors, so staff read settings quickly and confidently. Moreover, the SW30TG supports strict accuracy and efficient workflow. The built-in temperature data recorder includes a real-time clock, and it allows teams to program solution warming between 90°F and 110°F.

The self-tune feature maintains temperatures within $\pm 1^\circ\text{F}$, and the calibration reminder keeps facilities compliant.

FEATURES

- Single Compartment
- Comes with a Standard Glass Door
- Open, the adjustable design adds convenience and flexibility
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads

SPECIFICATIONS

- Model Number: SW30TG
- Part Number: 14BSW30243
- Dimensions: 30"W x 20-5/8"D x 60"H
- Shelves: 3
- Amps: 7.5

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online](#) →