

Blickman Solution/Blanket Warming Cabinet | 7921TS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7921TS solution/Blanket Warming Cabinet with Stainless Steel Door. The Blickman 7921TS Warming Cabinet delivers reliable heating performance for blankets or solutions. This single-door unit stands 74.5 inches tall, and it offers generous interior capacity for busy medical environments.

The open, adjustable shelf design increases flexibility, and it supports a wide range of warming needs. Meanwhile, the user can easily program solution-only warming between 90°F and 110°F for precise control.

Additionally, the Blickman 7921TS cabinet features simple, top-mounted controls that allow quick temperature adjustments. The PID digital controllers include a tactile stainless steel dome membrane and full numeric keypads, which enhance accuracy and usability.

FEATURES

- Single Compartment
- Single Stainless Steel Door
- 3 Shelves, open, adjustable design adds convenience and flexibility
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads

SPECIFICATIONS

- Model Number: 7921TS
- Part Number: 14B7921200
- Width: 30"
- Depth: 26-5/8"
- Height 74-1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)