

Amvex 2 Mode Vacuum/Suction Regulator | VR-C2U2



Stock Photo - Options May Vary

FACTORY NEW

Amvex 2 Mode Vacuum/Suction Regulator Series: VR-C2U2. Amvex® Vacuum Regulators provide the most advanced technology in the healthcare industry. The modularity allows for an easy upgrade from Analog to Digital right in the facility and from Continuous to Continuous Intermittent models.

Our units also allow for easy serviceability and maintenance. Amvex VR-C2U2 Series Features: Continuous 2 Mode (Off/Reg), 0-300 mmHg Gauge, Removable backplate and gasket for easy cleaning, Reinforced inlet and outlet ports to prevent breakage, Easily serviced and cleaned with cold sterilant, and Double O-Ring regulating module.

Your Choice of Outlet Connection (DISS, Ohmeda, PB, etc.), and comes with a 5 Year Limited Warranty.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)