

Blickman Solution/Blanket Warming Cabinet | 7924TG



Stock Photo - Options May Vary

FACTORY NEW

Blickman Solution/Blanket Warming Cabinet: 7924TG. The Blickman 7924TG Warming Cabinet stands 74.5 inches tall and delivers powerful dual-zone heating. It includes two separately heated compartments, so your team warms blankets and solutions at the same time. The key and passcode lockout system controls temperature settings and prevents unauthorized adjustments.

In addition, the cabinet uses a self-closing door with a magnetic gasket that improves efficiency and convenience. Furthermore, the 7924TG offers unrestricted temperature settings in both compartments. This flexibility allows your facility to manage diverse warming requirements with ease.

The PID digital controllers feature a tactile stainless steel dome membrane and full numeric keypads, which improve precision and user control.

FEATURES

- Double Compartments
- Comes with a Standard Double Glass Door
- Independently Controlled Compartments
- 3 Shelves, open, adjustable design adds convenience and flexibility
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)

SPECIFICATIONS

- Model Number 7924TG
- Part Number 14B7924243
- Width 30"
- Depth 36-5/8"
- Height 74.5"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com