

# Amvex Dual Oxygen Flowmeter with 1/8" FNPT | YB-FM-15UO-F2



Stock Photo - Options May Vary

## FACTORY NEW

Amvex Dual Oxygen Flowmeter with 1/8" FNPT: YB-FM-15UO-F2. Amvex® offers a complete line of Flowmeters from Neonatal to Adult. Flowmeters are available in the Standard Thorpe style (available in Oxygen, Medical Air, Heliox, CO<sub>2</sub>, and N<sub>2</sub>O) and Dial Style (available in Oxygen and Medical Air).

Products may appear slightly different from what is shown. The Amvex Flowmeter is a User-Friendly Device. Also Available in 0-3 LPM and 0-8 LPM.

Color-Coded Inner Tube and Labeling, No Strip "D" Shaped Knob, and is available in Oxygen, Medical Air, Heliox, CO<sub>2</sub>, and N<sub>2</sub>O. It also features back-pressure compensation to maintain flow accuracy despite restrictions and includes an inlet filter and extension.

## FEATURES

- Thorpe Style
- Amvex Dual Oxygen Flowmeter with 1/8th FNPT
- Model Number YB-FM-15UO-F2
- 0-15 LPM Dual Oxygen Flowmeter
- Available in 0-3 LPM and 0-8 LPM

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)