

# Precision Medical Chrome Flowmeter T Branch



Stock Photo - Options May Vary

## FACTORY NEW

The Precision Medical Chrome Flowmeter T Branch assembly expands medical gas delivery from a single wall outlet. It connects two chrome flowmeters for efficient clinical distribution. Built from lightweight aluminum, the assembly supports durability without added bulk.

Moreover, the T Branch delivers oxygen or medical air to two separate care areas with consistent flow. The optional power take-off can support advanced breathing equipment. Clinicians can connect ventilators or BiPAP machines directly to the demand valve.

This makes setup faster and reduces the number of extra connections. The centered design lines up neatly with the wall outlet, and right or left offset options are also available for custom setups.

## FEATURES

- Highly accurate needle valves ensure precise flow-settings
- Supply pressure of 50 psi
- Air and oxygen models available
- Green knob marked with "O2" for easy identification
- Durable polycarbonate housing

## SPECIFICATIONS

- Model: 1MFA1XXX
- Flow Range: 0-15 l/min
- Graduations: 0.5 l/min (0-5), 1.0 l/min (5-15)
- Accuracy
- 0-5 l/min  $\pm$  0.25 l/min
- 6-15 l/min  $\pm$  0.50 l/min
- Model: 1MFA11XX
- Flow Range: 2-26 l/min

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)