

AMVEX Oxygen Flowmeter with DISS Hand Tight Adapter | FM-15UO-DH

FACTORY NEW

AMVEX Oxygen Flowmeter with DISS Hand Tight Adapter: FM-15UO-DH. Amvex® offers a complete line of Flowmeters from Neonatal to Adult. Flowmeters come in Standard Thorpe style (Oxygen, Medical Air, Heliox, CO2, N2O) and Dial Style (Oxygen, Medical Air).

All flowmeters include a standard 1/8" FNPT connection, with options to fit different adapters and accessories for customization. MR conditional options are also available.

The Amvex Thorpe Flowmeters are CE-marked and designed to meet all current CGA, HCA, ISO, and FDA requirements. Amvex Flowmeters are available in USA and ISO colors.



Stock Photo - Options May Vary

FEATURES

- Thorpe style
- DISS Hand Tight Adapter
- Color-coded inner tube and labeling
- Easy to service
- No strip "D" shape knob

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)