

# AMVEX Oxygen Flowmeter with Ohmeda/Ohio Adapter | FM-15UO-OH



Stock Photo - Options May Vary

## FACTORY NEW

The Amvex FM-15UO-OH is a Thorpe-style oxygen flowmeter with Ohmeda/Ohio adapter connection, designed for reliable gas delivery in medical applications. It provides a flow range of 0-15 LPM with color-coded inner tube and labeling for easy identification. The flowmeter features backpressure compensation to maintain accurate flow rates despite downstream restrictions.

This model includes an inlet filter and inlet extension for enhanced protection and functionality. The no-strip "D" shaped knob design provides user-friendly operation and durability. The flowmeter connects via standard 1/8" FNPT connection and can be customized with different adapters as needed.

The Amvex flowmeter meets all current CGA, HCA, ISO, and FDA requirements and carries CE-marking certification. It is available in both USA and ISO color standards and is designed for easy servicing and maintenance throughout its operational life.

## FEATURES

- Thorpe style
- Color-coded inner tube and labeling
- Easy to service
- No strip "D" shape knob
- Includes an inlet filter and inlet extension
- Backpressure compensated so flow accuracy is not affected by restrictions

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)