

# Amvex Oxygen Flowmeter with Puritan Bennett Adapter | FM-15UO-PB

## FACTORY NEW

The Amvex Oxygen Flowmeter Model FM-15UO-PB is a Thorpe-style flowmeter with a 0-15 LPM oxygen delivery range and Puritan Bennett adapter. It features color-coded inner tube and labeling for easy identification, along with an inlet filter and inlet extension for enhanced functionality.



*Stock Photo - Options May Vary*

This flowmeter is backpressure compensated, ensuring flow accuracy is not affected by downstream restrictions. The design includes a standard 1/8" FNPT connection and is engineered for easy servicing with a no-strip "D" shape knob.

The unit is CE-marked and designed to meet current CGA, HCA, ISO, and FDA requirements. It is available in both USA and ISO color standards and can be configured with various adapters and accessories to meet specific clinical needs.

## FEATURES

- Thorpe style
- 0-15 LPM Oxygen Flowmeter
- Color-coded inner tube and labeling
- Easy to service
- No strip "D" shape knob
- Includes an inlet filter and inlet extension
- Backpressure compensated so flow accuracy is not affected by restrictions

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)