

Amvex Pediatric 2 Mode Vacuum/Suction Regulator | VR-P2U2



Stock Photo - Options May Vary

FACTORY NEW

The Amvex Pediatric 2 Mode Vacuum/Suction Regulator VR-P2U2 Series delivers advanced healthcare technology for safe, consistent clinical performance in pediatric care environments. The unit features a clear 0-150 mmHg gauge for instant pressure monitoring, enabling clinicians to achieve accuracy and reliability. Pediatric 2 Mode operation with Reg and Off settings provides controlled suction delivery.

One housing supports both Continuous and Intermittent regulation, simplifying inventory management and reducing facility costs. The removable backplate and gasket allow fast cleaning, while reinforced inlet and outlet ports resist breakage and extend equipment lifespan. A double O-ring regulating module improves sealing and performance.

The regulator supports easy servicing using cold sterilant for routine maintenance. Users can choose outlet connections including DISS, Ohmeda, or PB to ensure seamless integration across diverse medical gas systems. The unit includes a Limited 5-Year Warranty for long-term confidence.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)