

Amvex Vacuum / Suction Regulator - Cont/Inter | VR-CIU2



Stock Photo - Options May Vary

FACTORY NEW

The Amvex Vacuum/Suction Regulator VR-CIU2 Series combines digital display technology with the reliability of traditional analog gauges. The regulator features no moving parts, requires no calibration, and delivers $\pm 1\%$ full-scale accuracy for precise suction control.

The large digital display is readable from up to 30 feet away, with color-coded bars indicating Low, Medium, and High pressure ranges. The gauge operates across a 0-300 mmHg range and accommodates both continuous and intermittent suction applications within the same housing.

Constructed with durable and flexible polymer plastic casing, the unit includes a removable backplate and gasket for straightforward cleaning and maintenance. Multiple outlet connection options are available to match existing medical gas systems.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)