

# Amvex Vacuum/Suction Regulator Continuous High Vac | VR-CHU2



Stock Photo - Options May Vary

## FACTORY NEW

The Amvex VR-CHU2 Series is an analog vacuum/suction regulator designed for continuous high vacuum performance in healthcare settings. It features a 0-300 mmHg gauge and a modular design that allows on-site upgrades from analog to digital or from continuous to intermittent models, supporting long-term functionality while reducing maintenance costs.

The regulator includes a removable back plate and gasket for simplified cleaning and maintenance. Reinforced inlet and outlet ports enhance durability in high-demand environments, while a double O-ring regulating module ensures reliable performance. The device is compatible with cold sterilants for easy sterilization.

Multiple outlet connection options are available, including DISS, Ohmeda, and PB configurations to accommodate different facility setups. The product includes a 5-year limited warranty.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)