

Baxter 6201 FloGuard Infusion Pump - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Baxter 6201 FloGuard Infusion Pump is fully biomedically inspected and ready for immediate use in hospitals, clinics, and other healthcare facilities. It provides reliable performance for IV fluid management with Anti-Free Flow Protection for added safety during operation.

The pump offers two flexible flow rate settings: 1-99.9 mL/hr in 0.1 mL increments or 1-1,999 mL/hr in 1 mL increments. These options enable precise and accurate fluid delivery for varied patient treatment requirements.

Every unit includes a warranty and is compatible with universal sets. The Baxter 6201 FloGuard delivers dependable infusion support whether replacing older equipment or expanding facility inventory.

FEATURES

- Completely Biomedically Inspected and Patient Ready
- KVO 5 ML/HR or Lower Programmed Rate
- Anti-Free Flow Protection
- Universal Sets
- Warranty

SPECIFICATIONS

- Mounting: Pole
- Dimensions: 7.9 W x 11.4" H x 5.1" D
- Weight: 11.6 lbs
- Flow Rate 1-99.9 ML/HR in 0.1 ML /HR OR 1-1,999 ML/HR in 1 ML Increments
- KVO Rate: 5 ML/HR or Lower Programmed Rate

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)