

Baxter 6301 Dual Channel Infusion Pump - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Baxter 6301 is a dual-channel infusion pump designed for simultaneous delivery of fluids or medications in clinical settings. Each unit has been completely refurbished and biomedically inspected to meet high-quality standards. The pump supports accurate and safe delivery for adult and pediatric patients.

Key clinical features include dual volumetric IV delivery, Flow Check Occlusion Alarm, and Flow Rate Calculation for precision dosing. The pump offers Incremental Flow Rate control and Programmed Delivery Profiles to accommodate various infusion therapies. Additional capabilities include automatic restart once occlusion clears and automatic piggybacking for secondary medications.

The package includes a new battery, pole clamp, and user manual. The pump mounts to standard IV poles using the included clamp and is backed by warranty coverage.

FEATURES

- Dual Volumetric I.V. Pump
- Flow Check Occlusion Alarm
- Flow Rate Calculation
- Programmed Delivery Profiles
- Incremental Flow Rate

SPECIFICATIONS

- Dimensions: 13"W x 11.4" H x 5.1" D
- Weight: 17.9 lbs
- Rate: 0.1-1999 ml/hr
- Volume: 0.1-9999 ml
- Mounting: Pole (Pole Clamp is Included)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)