

Baxter InfusOR Syringe Infusion Device - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Baxter InfusOR Syringe Infusion Device is a refurbished pump used to give medicine through a vein during anesthesia. It works well for delivering painkillers, muscle relaxants, and drugs that affect blood pressure.

Its small size and reliability make it a great choice for outpatient surgeries where speed and accuracy are important. The InfusOR pump simplifies dose modifications without requiring multiple reprogramming steps.

The Baxter InfusOR allows clinicians to adjust settings easily at any point during the infusion process. The dial-based interface allows one-handed use, enabling fast adjustments when timing is critical during surgical procedures.

FEATURES

- For the administration of many intravenous agents given during anesthetic procedures
- Excellent device for outpatient surgical settings
- Modify dose without multiple reprogramming steps
- Adjusts settings easily at any time during infusion
- Dial settings can be changed with a single hand for an immediate response

SPECIFICATIONS

- Accuracy: Infusion linear rate $\pm 3\%$
- Batteries: Four "C" Alkaline Cells
- Typical Use Battery Life: 150 hours at low battery
- Operating Voltage: 5-7 volts
- Size: 9.2 x 4.5 x 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)