

# Biodex 240 MRI Stretcher – Completely Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Biodex 240 MRI Stretcher is a completely refurbished patient transport solution designed for MRI environments. Constructed from non-magnetic aluminum, it is MR Conditional for use in 3T MRI settings, ensuring safe compatibility during scanning. The lightweight design enables easy maneuverability while maintaining durability for reliable clinical use.

Equipped with 5-inch ball-bearing swivel casters and locking mechanisms, the stretcher provides controlled movement and secure positioning. Retractable guard rails offer patient protection, and the unit supports up to 500 pounds to accommodate a wide range of patient needs.

Each unit includes a new mattress and pad, comes fully refurbished and tested, and is backed by a 6-month warranty for long-term performance and reliability.

## FEATURES

- Heavy Duty MRI Stretcher
- Light Weight All Aluminum-Non-magnetic
- Perfect for Patient Transport
- Retractable Side Rails
- Easy to roll on 5" locking, swivel casters
- Includes New Pad
- Warranty

## SPECIFICATIONS

- Tabletop Dimensions: 72" L x 24" W
- Overall Dimensions: 74.5" L x 27.75" W x 34" H
- Casters: 5" locking, swivel
- Patient Capacity: 500 lb
- Compliance: MR Conditional
- Weight: 85 lbs
- Warranty: 6 months
- Mattress: New

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)