

Hausted MBC Mammography Biopsy Chair - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hausted MBC Mammography Biopsy Chair is a versatile unit designed for breast imaging and minor procedures. Its compact base fits easily into tight rooms with bulky imaging equipment, and the backrest adjusts to a full 90° upright position to streamline mammography exams. The chair smoothly transitions patients from sitting to supine positions while maintaining comfort and stability.

Foot pedals on both sides allow effortless height and Trendelenburg adjustments without manual intervention. Side rails fold down easily using quick-release levers, providing better access to lesion sites during biopsies or treatments. The built-in brake and steering system features easy-to-reach levers on the chair base for convenient operation.

The stainless steel base cover is easy to wipe clean for hygiene maintenance. An adjustable, fold-down footrest enhances lower body support and stability during longer procedures. The chair includes a built-in chart holder and push bar to facilitate movement and document organization.

FEATURES

- The unit's height and Trendelenburg positions can be effortlessly adjusted using foot pedals located on both sides of the base, allowing for hands-free operation.
- Side rails fold away easily via conveniently placed release levers, ensuring optimal access to the lesion site.
- A brake and steering system is controlled through side-mounted levers on the chair base.
- The seamless stainless steel base cover is designed for easy cleaning.
- Includes an integrated chart holder and a push bar for added convenience.
- An adjustable, fold-down footrest enhances lower body comfort and stability.
- Includes New Pad Set

SPECIFICATIONS

- Height Range: High: 31.5" (80 cm), Low: 23.5" (60 cm).
- Standard Trendelenburg: 0 - 10°.
- Patient Surface Width: 22" (56 cm).
- Overall Length: Supine: 75" (191 cm), Upright: 42.5" (108 cm).
- Wheel Base: 25" x 22" (63 cm x 55 cm).
- Pneumatic-Assist Backrest: 0 - 90°.
- Side Rail Height: 10.75" (27 cm).
- Maximum Weight Capacity: 325 lbs (evenly distributed)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com