

Stryker 1080 Full Length C-Arm Access Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Refurbished Stryker 1080 Stretcher features a U-shaped base that offers full-length C-arm access, improving imaging capabilities and allowing full access to the patient during procedures. Fold-down side rails stay securely out of the way when needed, helping caregivers work more efficiently. Dual pneumatic backrests on both ends allow positioning from either side, reducing the need to move patients during procedures.

The stretcher includes a radiolucent top for optimal imaging and supports Trendelenburg and Reverse Trendelenburg positioning. Staff can safely lower or raise the patient using dual-end hydraulic foot controls, while a four-wheel steel ring brake system with dual-end activators helps secure the treatment platform.

With a 500 lb patient weight capacity and manual height adjustment ranging from 24 to 34 inches, the Stryker 1080 is designed for versatile clinical use. The unit comes with a new 2-inch pad and includes four IV receptacles for patient support.

FEATURES

- Radiolucent Litter Top
- Dual Fiberesin Pneumatic Backrest
- Dual-End Hydraulic Foot Controls
- Fold-Down Side-Rails
- Four-Wheel, Steel Ring Brake System with Dual-End Activators
- Four IV Receptacles
- 2" Thick Stretcher Pad
- 6" Non-Conductive Casters

SPECIFICATIONS

- Weight Capacity: 500 lbs
- Overall Bed Length: 87"
- Overall Bed Width: 30.5"
- Manual Height Adjustment: 24" - 34"
- Patient Surface Dimensions: 25.5" x 84"
- Side-Rail Dimensions: 13" x 55"
- Warranty: 6 Months Parts
- Completely Patient Ready

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)