

# Mindray Datascope Passport V Patient Monitor- Refurbished



*Stock Photo - Options May Vary*

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray Datascope Passport V Patient Monitor – Refurbished. The Mindray Datascope Passport V is a great fit for most departments in an ASC or hospital. Additionally, it features a 12.1-inch color LCD screen and battery backup for continuous monitoring.

Furthermore, this multi-parameter monitor includes ECG, NIBP, RESP, SpO2, Temp, and an internal thermal recorder for essential patient data. Moreover, this monitor offers 96 hours of graphic and list trends, along with visual and audible alarms. As a result, staff can review patient history and stay alert to changes.

Specifically, drug calculations and remote viewing with multi-bed alarm notifications add further support.

## FEATURES

- 12.1" Color LCD
- ECG 3/5 Lead
- Non-Invasive Blood Pressure
- Impedance Respiration
- Masimo SET SpO2

## SPECIFICATIONS

- Display Size (Screen): 12.1"
- Battery: Rechargeable, Lithium-Ion
- 96 Hours of Graphic and List Trends
- 128 Minutes of Oxy-CRG Trends
- 3-Trace Recorder

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)