

Mindray Datascope Passport V Patient Monitor- Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray Datascope Passport V Patient Monitor – Refurbished.

The Mindray Datascope Passport V is perfect for most departments in an ASC or hospital.

Furthermore, it features a 12.1" color LCD screen and a battery backup for continuous monitoring. Additionally, the multi-parameter monitor includes ECG, NIBP, RESP, SpO2, Temp, and an internal thermal recorder, providing essential data for patient care.

Additionally, the Mindray Datascope Passport V offers 96 hours of graphic and list trends, visual and audible alarms, drug calculations, and remote viewing with multi-bed alarm notifications. This Passport V also supports device connectivity via DIAP and Ethernet, and includes a rechargeable lithium-ion battery for extended use.

FEATURES

- 12.1" Color LCD
- ECG 3/5 Lead
- Non-Invasive Blood Pressure
- Impedance Respiration
- Masimo SET SpO2

SPECIFICATIONS

- Display Size (Screen): 12.1"
- Battery: Rechargeable, Lithium-Ion
- 96 Hours of Graphic and List Trends
- 128 Minutes of Oxy-CRG Trends
- 3-Trace Recorder

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)