


Edan iM70 Patient Monitor – New



Stock Photo – Options May Vary

FACTORY NEW

Edan iM70 Patient Monitor – New. The Edan iM70 Patient Monitor offers a wide range of parameters, from invasive blood pressure to cardiac output. Additionally, it also covers capnography and anesthesia gases.

Furthermore, this versatility helps meet many clinical needs across different settings, including general wards, cardiac units, ICUs, and operating theaters. Moreover, the iM70 comes with standard parameters such as 3- or 5-lead ECG, RESP, and NIBP. Additionally, it includes Edan SpO2, dual TEMP, PR, and HR.

As a result, staff get a full picture of patient status in one place. Specifically, give us a call to discuss other available options for this monitor.

FEATURES

- 12" Color TFT LCD Display
- 120 Hour Trend Storage
- ECG 3 or 5 Lead Available
- Non-Invasive Blood Pressure
- SpO2

SPECIFICATIONS

- Display Size (Screen): 12"
- Dimensions: 12.9" L x 6.2" D x 11.2" H
- Weight: 12.1 lbs.
- Battery: Rechargeable, Lithium-Ion
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)