

Edan iM8 Patient Monitor with CO2 – New



Stock Photo – Options May Vary

FACTORY NEW

Edan iM8 Patient Monitor with CO2 – New – The Edan iM8 puts reliable monitoring technologies at your fingertips. With its accurate readings and optimized performance, the iM8 can meet basic monitoring needs, bringing you easy to high-quality healthcare for all patients. The Edan G2 CO2 (sidestream) – Superior water trap design for accurate monitoring iCARB™ algorithm with Intelligent CO2 pseudo wave identification technology In Addition, the Edan iM8 Features: 12" LCD Screen Pacemaker Detection HL7/LAN Connection to EMR/Central Station Built-in Rechargeable Lithium Battery Pitch Tone (pulse-tone modulation) Dual-Mode Anti-Interference Pulse Oximetry Defibrillation and Electrosurgical Interference Protection Unique iSEAP algorithm optimized for arrhythmia patients Are you looking for additional monitors and accessories? Click the following link to view our complete catalog of PATIENT MONITORS FOR SALE | NEW, USED & REFURBISHED.

FEATURES

- 12.1" Color TFT LCD Display
- Standard Edan G2 CO2 (sidestream) or optional - Respiration CO2 (mainstream/sidestream)
- 96 Hour Trend Storage
- ECG 3 or 5 Lead Available
- Non-Invasive Blood Pressure

SPECIFICATIONS

- Display Size (Screen): 12.1"
- Dimensions: 12.6" L x 5.9" D x 10.4" H
- Weight: 11 lbs.
- Battery: Rechargeable, Lithium-Ion
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)