

# Edan iM60 Patient Monitor with CO2 - New



Stock Photo - Options May Vary

## FACTORY NEW

Edan iM60 Patient Monitor with CO2 - New - The iM series offers a comprehensive array of parameters, spanning from invasive blood pressure to cardiac output, and from capnography to anesthesia gas.

This versatility allows it to meet the diverse clinical needs across various settings.

## FEATURES

- 10.4" Color TFT LCD Display
- Standard Edan G2 CO2 (sidestream) or optional - Respironics CO2 (mainstream/sidestream)
- 120 Hour Trend Storage
- ECG 3 or 5 Lead Available
- Non-Invasive Blood Pressure

## SPECIFICATIONS

- Display Size (Screen): 10.4"
- Dimensions: 11.9" L x 6.3" D x 10" H
- Weight: 9.9 lbs.
- Battery: Rechargeable, Lithium-Ion
- Power Requirement: 120v

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)