

Mindray DPM 6 Patient Monitor with CO2 – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Mindray DPM 6 Patient Monitor w/ CO2 is designed to meet the needs of the ICU and operating room, offering advanced features such as multi-vector Arrhythmia and ST analysis, Cardiac Output (CO), Continuous Cardiac Output (CCO/SVO2), CO2 monitoring, up to 8 invasive blood pressure measurements, and 5-agent multi-gas analysis with automatic identification and data storage capabilities. The Mindray DPM 6 features comprehensive monitoring with advanced data management and an intuitive, customizable interface.

It integrates seamlessly with hospital networks, enabling easy access to patient data. The high-resolution display clearly shows vital parameters.

Its robust alarm system identifies critical changes, allowing for timely interventions and improved patient care.

FEATURES

- Robust measurement set for mid to high-acuity environments
- Enhanced data storage and recall capabilities
- 3, 5, and 12-lead ECG
- Arrhythmia and ST analysis (optional)
- Up to 8 invasive blood pressures (optional)

SPECIFICATIONS

- Size: 11.7" (H) x 13.2" (W) x 7.3" (D) 29.7 cm (H) x 33.6 cm (W) x 18.6 cm (D)
- Weight: 15.8 lbs (7.2 kg)
- Including 12.1" touchscreen display, MPM, and ECG, NIBP, SpO2 accessories
- AC Voltage: 100 to 240VAC, 50/60Hz
- Current: 2.5 to 1.4A

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)