

Mindray DPM 6 Anesthesia Monitor with AG Module - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray DPM 6 Anesthesia Monitor with AG Module - Refurbished. The Mindray DPM 6 Patient Monitor is well-suited for most departments in an ASC or hospital.

Additionally, the DPM 6 features a 12.1" Color Touchscreen Display and Battery Backup. Additionally, this Multi-Parameter Monitor comes standard with ECG, NIBP, RESP, SpO2, Temp, and Internal Thermal Recorder.

Furthermore, this Mindray DPM 6 includes an anesthesia gas module (AG) for monitoring 5 agent anesthesia gases and end tidal CO2.

FEATURES

- 12.1" Color Touchscreen
- ECG 3/5 (12 lead optional)
- Non-Invasive Blood Pressure
- SpO2
- Continuous Temperature

SPECIFICATIONS

- Display Size (Screen): 12.1"
- Battery: Rechargeable, Lithium-Ion (2)
- 120 Hours of Trend data at 1 min resolution/ 1 Hour at 1 sec resolution
- 3-Trace Thermal Array Recorder
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)