

Mindray DPM 7 Patient Monitor with AG Module - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray DPM 7 Patient Monitor with AG Module is a powerful and adaptable monitoring system, ideal for mid to high acuity settings such as the ICU and OR. It offers advanced capabilities including multi-vector arrhythmia and ST analysis, cardiac output monitoring, and up to eight invasive blood pressure channels.

The Anesthesia Gas Module provides 5-agent multi-gas analysis with automatic agent identification and data storage, supporting accurate anesthetic monitoring. The large 17" high-resolution touchscreen enhances visibility and usability, while the built-in 3-trace recorder simplifies documentation.

Designed with flexibility in mind, the Mindray DPM7 offers independent display configurations for customized patient data views, making it ideal for anesthesiologists, surgeons, and perfusionists.

FEATURES

- Robust measurement set for mid to high acuity environments
- Enhanced data storage and recall capabilities
- 3/5 Lead
- Arrhythmia and ST analysis (optional)
- Up to 8 invasive blood pressures (optional)

SPECIFICATIONS

- Dimensions: 14.6" (370 mm) H x 15.7" (400 mm) W x 7.6" (193 mm) D
- Weight: 31.9 lbs. (14.5 kg) including display, MPM, AG module, 2 lithium batteries, recorder, CF components
- 12.1" touchscreen display, MPM, and ECG, NIBP, SpO2 accessories
- AC Voltage: 100 to 240VAC, 50/60Hz
- Current: 2.5 to 1.4A

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)