

Mindray DPM7 Patient Monitor - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Mindray DPM7 is a multi-parameter patient monitor designed for ICU and OR environments. It features a 17-inch high-resolution touchscreen with a unique split-screen capability that allows simultaneous viewing of real-time data and 5-lead ECG. The monitor supports comprehensive patient assessment including arrhythmia detection, 3/5/12-lead ECGs, respiration, non-invasive blood pressure, ST analysis, cardiac output, and mixed venous oxygen saturation.

The DPM7 includes auto-ID functionality and data storage features to enhance workflow and record-keeping. Optional modules for up to 8 invasive pressures, CO2 monitoring, and five-agent gas analysis are available separately. The unit comes equipped with new accessories including ECG patient cable and lead wires, SpO2 finger sensor, NIBP cuff and hose, lithium-ion battery, power cord, printer paper, and user manual.

FEATURES

- Clear 17" high-resolution color touchscreen
- Comes standard with lithium-ion batteries
- Customizable touchscreen quick-action keys for efficiency
- Large numerics display mode for increased visibility

SPECIFICATIONS

- Physical Dimensions: 14.6" (370mm) H x 15.7" (400mm) W x 7.6" (193cm) D
- Weight: > 31.9lbs (14.5 kg) including display, MPM, AG module, 2 lithium batteries, recorder, CF components
- Display: 17" (432 mm)
- Input Voltage: 100 to 240VAC
- Operating Temperature: 41° - 104°F (5° to 40°C)
- Operating Humidity: 15% to 95%, non-condensing
- Network: 2 RJ45 network connectors, 100 BASE-TX
- Power Supply: 1 AC power connector

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)