

Mindray Datascope Passport V Anesthesia Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray Datascope Passport V Anesthesia Monitor – Refurbished. This monitor fits well in most surgery departments of an ASC or hospital. The Passport V with Gas Module 3 covers all standard Passport V parameters.

In addition, it can track CO2 and five anesthesia agents during surgery. This refurbished Passport V has a 12.1-inch color LCD screen and battery backup. The multi-parameter monitor comes standard with ECG, NIBP, RESP, SpO2, temperature, and an internal thermal recorder.

Additionally, the included Gas Module 3 tracks end-tidal CO2 and five anesthetic agents with automatic ID.

FEATURES

- 12.1" Color LCD
- ECG 3/5 Lead
- Impedance Respiration
- Masimo SET SpO2
- Continuous Temperature

SPECIFICATIONS

- Display Size (Screen): 12.1"
- Battery: Rechargeable, Lithium-Ion
- 96 Hours of Graphic and List Trends
- 128 Minutes of Oxy-CRG Trends
- 3-Trace Recorder

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)