

Mindray Datascope Passport V Anesthesia Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Mindray Datascope Passport V is a full-featured anesthesia monitor designed for surgery departments in hospitals and ambulatory surgical centers. It includes a 12.1" color LCD display and comes standard with ECG, NIBP, RESP, SpO2, temperature monitoring, and an internal thermal recorder. The integrated Gas Module 3 enables monitoring of end-tidal CO2 and five anesthetic agents with automatic identification.

This multi-parameter monitor provides 96 hours of graphic and list trends, visual and audible alarms, and drug calculation capabilities. It supports up to 8 waveforms with automatic configuration and dual invasive pressure monitoring for comprehensive patient assessment during procedures.

The unit features a rechargeable lithium-ion battery for portability and includes device connectivity via DIAP and Ethernet for remote viewing and multi-bed alarm notifications. A 3-trace recorder and extensive trend storage make it suitable for detailed perioperative monitoring and documentation.

FEATURES

- 12.1" Color LCD
- ECG 3/5 Lead
- Masimo SET SpO2
- Lithium-Ion Battery Support
- Device Connectivity via DIAP and Ethernet

SPECIFICATIONS

- Display Size (Screen): 12.1"
- Battery: Rechargeable, Lithium-Ion
- 96 Hours of Graphic and List Trends
- 3-Trace Recorder
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)