

Datascope Passport 2 Anesthesia Monitor with Gas Module – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Datascope Passport 2 Anesthesia Monitor with Gas Module: Refurbished. The Datascope Passport 2 Anesthesia Monitor allows healthcare professionals to quickly assess patient status at a glance of the screen. Its user-friendly design and clear display make it an essential tool for effective patient monitoring during procedures.

Key Features and Benefits of the Passport 2: User-friendly design and menus: Simplifies operation and enhances efficiency for healthcare providers. Large 10.4" bright color display: Ensures visibility, even from a distance, for easy monitoring in busy environments.

Monitoring of ECG, NIBP, SpO2, and Temp: Provides continuous tracking of essential patient vitals, ensuring comprehensive care.

FEATURES

- Multi-Parameter
- Measures ECG, NIBP, SpO2, Temp
- Gas Module for Co2 and AA
- Auto Agent Identification
- 3 or 5 Lead ECG

SPECIFICATIONS

- Dimensions: ~14"H x 12"W x 11"D
- Weight: ~24 Lbs
- Display Size (Screen): 10.4"
- Recorder Paper: 80 mm, Thermal
- Battery: Rechargeable, Sealed Lead

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)