


# Mindray Datascope Passport 2 Anesthesia Monitor with Gas Module – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Datascope Passport 2 Anesthesia Monitor with Gas Module – Refurbished. This monitor lets staff check patient status at a glance.

Its easy-to-use design and clear display make it a key tool for monitoring during procedures. The large 10.4-inch bright color display stays visible even from a distance, which helps in busy settings.

In addition, it tracks ECG, NIBP, SpO2, and temperature for full, ongoing care. The gas module attaches under the monitor to measure EtCO2, FiO2, O2, N2O, and five anesthetic agents during procedures.

## FEATURES

- Multi-Parameter
- Measures ECG, NIBP, SpO2, Temp
- Gas Module for Co2 and AA
- Auto Agent Identification
- 3 or 5 Lead ECG

## SPECIFICATIONS

- Dimensions: ~14"H x 12"W x 11"D
- Weight: ~24 Lbs
- Display Size (Screen): 10.4"
- Recorder Paper: 80 mm, Thermal
- Battery: Rechargeable, Sealed Lead

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)