

Datascope Trio Patient Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Datascope Trio Patient Monitor – Refurbished. The Datascope Trio Patient Monitor – Refurbished is a versatile, cost-effective solution for patient monitoring. Its compact design makes it ideal for both bedside and transport use. Healthcare professionals can rely on its real-time monitoring of essential vital signs in various clinical settings, from hospitals to emergency response situations.

It provides an efficient way to monitor patients in fast-paced, demanding environments. Equipped with an 8.4" color display, the monitor offers clear and accurate readings. Healthcare providers can easily assess patient conditions without unnecessary complexity.

It includes all necessary features for basic monitoring, including ECG (Electrocardiogram), SpO2 (Oxygen Saturation), NIBP (non-invasive blood pressure), impedance respiration, and temperature readings.

FEATURES

- 8.4" TFT Color Display
- ECG 3/5 lead (4 traces with cascade)
- Impedance Respiration
- Non-Invasive
- Blood Pressure

SPECIFICATIONS

- Dimensions: Approximately 9"H x 9.5"W x 7"D
- Weight: Approximately 11 lbs
- Display Size (Screen): 8.4"
- Recorder Paper: 50mm, Thermal
- Recorder Traces: 2

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)