

GE Dash 2000 Patient Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

GE Dash 2000 Patient Monitor – Refurbished. The refurbished GE Dash 2000 Patient Monitor is a compact, affordable, and dependable device for continuous patient monitoring. Additionally, this monitor delivers fast, accurate data to support quality patient care.

Furthermore, it works well in both bedside and transport settings. Moreover, this monitor includes the GE EK-Pro arrhythmia detection program and GE Dinamap SuperSTAT non-invasive blood pressure. As a result, staff get reliable readings for quick patient checks.

Specifically, staff can choose Masimo or Nellcor SpO2 technology for accurate oxygen saturation readings. In addition, the Alaris Turbo Temp thermometer adds quick, reliable temperature monitoring.

FEATURES

- Functions: ECG, NIBP, PR, SpO2 plus Temp
- Easy to Use Display
- Compact Sizes
- Rechargeable Battery
- Visual and Audible Alarms

SPECIFICATIONS

- Display: 5.8" LCD
- Dimensions: 8.5" H x 10.2" W x 7.9" D
- Weight: 11.5 lbs
- Battery Run Time: Approx. 3.5 hours
- Charge Time: 3 Hours Max

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)