

GE Dash 2500 Patient Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

GE Dash 2500 Patient Monitor – Refurbished. The GE Dash 2500 Patient Monitor is a reliable, affordable device built for speed and accuracy. Additionally, this monitor helps staff manage patient care with precision.

Furthermore, the GE EK-Pro arrhythmia program helps staff spot heart rhythm issues quickly and clearly. Moreover, the GE Dinamap SuperSTAT gives fast, non-invasive blood pressure readings for quick patient checks. As a result, staff get reliable results without added wait time.

Specifically, staff can choose Masimo or Nellcor SpO2 for precise oxygen saturation readings. In addition, the Alaris Turbo Temp gives accurate, efficient temperature readings for full patient monitoring.

FEATURES

- Functions: ECG, NIBP, PR, SpO2 plus Temp
- Easy to Use Display
- Compact Size
- Rechargeable Battery
- Navigator Knob for Easy Configurations

SPECIFICATIONS

- Dimensions: 8.7"H x 14.1"W x 6.7"D
- Weight: 12 lbs
- Battery Run Time: Approx. 3 Hours
- Charge Time: 4 Hours Max.
- Power Requirement: 120V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)