

GE Dash 3000 Patient Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

GE Dash 3000 Patient Monitor – Refurbished. The GE Dash 3000 Patient Monitor is a reliable, affordable device built for accurate, efficient patient care. Additionally, the GE EK-Pro arrhythmia program helps staff spot heart rhythm issues quickly and clearly.

Furthermore, the GE Dinamap SuperSTAT gives fast, accurate blood pressure readings for quick checks. Also, this combination of tools supports faster, more confident clinical decisions. Moreover, this monitor offers Masimo or Nellcor SpO2 for precise oxygen saturation monitoring.

As a result, staff get reliable readings across many patient types. Specifically, the Alaris Turbo Temp gives accurate, efficient temperature readings. In addition, this monitor helps healthcare providers manage patient care effectively.

FEATURES

- Functions: ECG, NIBP, PR, SpO2 plus Temp
- Easy to Use Display
- Compact Size
- Rechargeable Battery
- Visual and Audible Alarms

SPECIFICATIONS

- Display: 8.4" LCD
- Dimensions: 10.25"H x 11" W x 8" D
- Weight: 11.2 lbs
- Recorder Paper: 50mm, Thermal
- Battery Run Time: Approx. 4 -5 Hours

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)