

# Birkova Lift Assist Leg Positioning System | LS-3500



*Stock Photo - Options May Vary*

## FACTORY NEW

The Birkova LS-3500 Lift Assist Leg Positioning System provides superior support and positioning for patients during medical procedures. The system features pneumatic lift assist with 120° rotation capability, handle-controlled abduction and adduction, and a twist handle for precise leg positioning that mimics natural hip motion. The complete system includes one set of pneumatic lift assist leg holders, one set of pressure management boot pads, and one set of siderail clamps.

## FEATURES

- Pneumatic Lift Assist
- 120° Rotation
- Handle Controlled Abduction and Adduction
- Mimics Natural Hip Motion
- Twist the Handle to Position the Leg
- 500 lbs. Weight Capacity
- Gel boot pads are available (see Options)

## SPECIFICATIONS

- Rotation: 120°
- Weight Capacity: 500 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)