

Birkova Knee Crutch Supports | LS-1700



Stock Photo - Options May Vary

FACTORY NEW

Birkova Knee Crutch Supports LS-1700.

Birkova Products LLC has been manufacturing and supplying a vast range of quality products for the medical profession since 1957, generating worldwide synergy in the Medical Industry with its products and services.

Birkova's LS-1700 Features: The LS-1700 supports are Height Adjustable with a 360° Full Rotation. Stainless Steel Support Bars and Attaches to Surgical Table using Super Clamps or Clamp-Tite™ Clark Socket (Sold Separately). Comes complete with Pressure management Pads and Restraint Straps.

FEATURES

- Height adjustable - 360° full rotation
- Stainless steel support bars
- Includes pressure management pads and restraint straps
- Attaches to the surgical table using Super Clamps or Clamp-Tite Clark Sockets (sold separately)

SPECIFICATIONS

- 360° Full Rotation
- Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)