

Amsco 2080 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Amsco 2080 Standard Series 3-piece replacement pad set is designed for models 2080M, 2080L, 2080RC, and 1080 surgical tables. This set is ideal for low-risk patients and short-duration procedures, providing reliable performance in any surgical environment.

Each pad features 2 inches of medium-firm high-density foam construction with an electrically conductive Lectrolite vinyl cover. The pads are fluid and stain resistant, X-ray permeable, and include anti-microbial properties for enhanced safety and durability.

The set includes a hook and loop kit for secure attachment and is latex-free for patient safety. These replacement pads offer comfort and reliable support during standard surgical procedures.

FEATURES

- For Amsco Table Models: 2080M, 2080L, 2080RC, and 1080
- For Low-Risk Patients
- For Short Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Hook & Loop Kit Included
- Anti-Microbial

SPECIFICATIONS

- Construction: Medium Firm High-density Foam
- Exterior Cover: Electrically Conductive Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)