

Amsco Steris 3080/3085 Bariatric Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Amsco Steris 3080/3085 Bariatric Series Replacement Pad Set is a 3-piece pad system designed specifically for bariatric and high-risk patients. It is compatible with Amsco Steris surgical table models 3080 RC, 3080 RL, 3080 SP, and 3085 SP.

The pad set features a 4-inch thick construction with a high-density firm foam core and an anti-shear memory foam layer. The cover is made of 3-way stretch, ultra-soft Fusion Fabric vinyl that resists fluid and stains.

This design reduces pressure points and skin shear during prolonged surgical procedures, providing maximum pressure relief and support for low to high-risk procedures involving obese or bariatric patients.

FEATURES

- Designed for the Comfort of Bariatric Patients
- For Low-High Risk Obese Patients
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free
- Hook & Loop Kit Included

SPECIFICATIONS

- Inside Layer: High-Density Comfort Foam
- Second Layer: Medium Firm High-Density Foam
- Third Layer: Prot-Tek memory foam
- Cover: 4-Way Stretch Fusion Fabric

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)